

Sunday, September 28, 2025

Designing oligonucleotides to disrupt lacZ gene in E coli DH5 alpha

Ali M.

---

Goal: Knock-out the lacZ gene.

Steps:

**Order two oligonucleotides** (instructions will be below).

**Anneal them** → mix forward and reverse oligos, heat and cool so they form a double-stranded insert.

**Clone into Cas9 plasmid** using Golden Gate (BbsI) or similar cloning strategy.

**Transform the plasmid into E. coli DH5α.**

Inside the cells, the plasmid will express **Cas9 + our gRNA**, which together target lacZ and cut it at the chosen site.

The cut will usually inactivate lacZ via small indels.

### Oligonucleotide design:

**Step 1:** Obtain the genbank/FASTA link for the E. Coli strand we're targeting. [This is the link.](#)

**Step 2:** Control-F to find the lacZ gene's coordinates. Note the CDS, or coding DNA sequence, maps to β-galactosidase, which is what we want to make non-functional.

```
gene      complement(1239256..1242168)
          /gene="lacZ"
          /locus_tag="NEB5A_06270"
CDS       complement(1239256..1242168)
          /gene="lacZ"
          /locus_tag="NEB5A_06270"
          /inference="EXISTENCE: similar to AA
sequence:RefSeq:WP_011378569.1"
          /note="forms a homotetramer; hydrolyzes lactose
disaccharide to galactose and glucose; converts lactose to
allolactose which is the natural inducer of the lac
operon; Derived by automated computational analysis using
gene prediction method: Protein Homology."
          /codon_start=1
          /transl_table=11
          /product="beta-D-galactosidase"
          /protein_id="A0072618.1"
```

**Step 3:** Using a python script, I created a new FASTA file that only contained nucleotides in the range of 1239256..1242168.

**Step 4:** Go to CHOPCHOP, and find the best protospacer. Here's a link to the CHOPCHOP output with the new FASTA file from Step 3.

**Step 5:** According to CHOPCHOP, our best protospacer is 'CTGACAATGGCAGATCCCAG'. So our forward oligo will be:

**CACC CTGACAATGGCAGATCCCAG**

And our reverse oligo will be:

**AAAC CTGGGATCTGCCATTGTCAG**

after we add the overhangs.

**Step 6:** Order the oligonucleotides on idtdna's website.

# 1 oligo\_fwd

ACTIONS

qty1GO

\$10.80

Product

25 nmole DNA Oligo

Expected Ship Date

10/1/2025

Purification

Standard Desalting

Guaranteed Yield

2.7 ODs =  
12 nmol =  
87.5 µgrams

Length

24

Sequence

CAC CCT GAC AAT GGC AGA TCC CAG

# 2 oligo\_rev

ACTIONS

qty1GO

\$10.80

Product

25 nmole DNA Oligo

Expected Ship Date

10/1/2025

Purification

Standard Desalting

Guaranteed Yield

2.8 ODs =  
12 nmol =  
88.2 µgrams

Length

24

Sequence

AAA CCT GGG ATC TGC CAT TGT CAG